

Liasys Clinical Chemistry System

Simple and Proven Design • Cost Effective • Easy to Use

**UP TO 330 TESTS
PER HOUR**

Technical Characteristics

LOADING SYSTEM:

REAGENTS RACKS:

4 Racks hold up to 33 reagents
2 Cooled racks up to 20 reagents
2 Racks with 13 reagent positions
and 14 Standards and Controls
2 Additional positions for diluents

SAMPLES RACKS:

5 Racks for samples loading
64 Positions for routine samples
14 Positions for STAT samples

FEATURES:




INTEGRATED ROBOTIC ARM

Air bath incubation system
60 WASHABLE CUVETTE WITH
QC MONITORING PROGRAM
USER CHANGEABLE CUVETTES
Capacitive liquid level sensing
High precision ceramic diluter
Side mounted ISE module
US Ready - ASTM Standards
Warranty Extension Plan

OPERATION:

WALK-AWAY:

Fully automated random access
Up to 330 Tests/Hour Bench-Top
GRAPHICAL INTERFACE
On-Board, 5-step cuvette
washing system
CONTINUOUS SAMPLE LOADING
14 STAT Positions
Optional Bar Code Readers
Users Help Menu



Best Value for your Laboratory

QC Monitored
On-Board Cuvette
Washing System

Durable Plastic
Cuvettes:
No-Breakages

Integrated
Robotic
Sampling Arm



TEST MENU

TEST PANELS

Basic Metabolic Panel
Comprehensive Metabolic Panel
Electrolyte Panel
Hepatic Panel
Lipid Panel
Renal Panel
User-defined panels and
ratio calculations available

INDIVIDUAL TESTS

General Chemistry Assays

Albumin
CO2
Direct Bilirubin
Total Bilirubin
Calcium
Creatinine
Glucose
Phosphorus
Magnesium
Total Protein
Urea Nitrogen (BUN)
Uric Acid

Electrolytes

Sodium
Potassium
Chloride

Enzyme Assays

ALT
AKL
Amylase
AST
CK(NAC)
Gamma GT
LDH

Lipid Assays

Cholesterol
Direct HDL
Direct LDL
Triglycerides

Specialty Assays*

Fructosamine
Hemoglobin Alc
hs-CRPWR
Homocysteine
Iron
Lipase
Microalbumin
Microprotein
Rheumatoid Factor

*Call Jant Pharmacal Corporation at
800.676.5565 for list of latest assays.



JANT PHARMACAL CORPORATION

16255 Ventura Blvd., #505, Encino, CA 91436

800.676.5565 818.986.8530 Fax 818.986.0235

www.accutest.net info@accutest.net

Accutest is a registered trademark of Jant Pharmacal Corporation