




accutest® Test Before You Treat!

The Accutest® Integrated Strep A Test provides results in just 5 minutes



CLIA Waived!

Accurate Results

- High clinical sensitivity at low colonization levels
- Built in procedural control
- External positive and negative quality controls provided

Simple Procedure

- Easy to perform and interpret immediately

Diagnosis

- Results in 5 minutes
- Facilitates patient treatment decisions in the same office visit
- Reduces risk of Strep A complications and spread of infection



Convenient

- Complete kit, no additional components required
- Room temperature storage
- Catalog No. ID440

The Integrated Strep A Rapid Test Device

Simple Procedure

- Add reagents
- Agitate swab, wait 1 minute
- Twist open valve
- Read results



JANT PHARMACAL CORPORATION

16255 Ventura Blvd., #505, Encino, CA 91436
800.676.5565 818.986.8530 Fax 818.986.0235
www.accutest.net info@accutest.net