

accutest[®] *Urine Analyzer* **500** *Your Choice for Value*



- High throughput: 514 strips/hour
- High luminosity cold source test system with 4 wavelengths provides superior sensitivity, specificity and accuracy
- Designed to work with Accutest URS 10 Urine Reagent Strips
- Stores up to 2000 test results
- Large LCD touch screen display
- Built-in high speed, low noise thermal printer
- Computer interface: RS-232 port
- Parallel printer interface
- External barcode reader may be connected to unit
- Includes waste collection system
- Lightweight: Less than 15 pounds
- Small footprint: 15"x14.9"x10.8"
- Catalog No. UA 750



JANT PHARMACAL CORPORATION

16255 Ventura Blvd., #505, Encino, CA 91436

800.676.5565 818.986.8530 Fax 818.986.0235

www.accutest.net info@accutest.net

Accutest is a registered trademark of Jant Pharmacal Corporation